



How I do it: en-bloc subaxial cervical laminectomy using a high-speed drill with a footplate attachment

Kyle B. Mueller¹ · Kenneth P. Mullinix² · Hector F. Bermudez³

Received: 18 November 2019 / Accepted: 25 November 2019 / Published online: 10 December 2019
© Springer-Verlag GmbH Austria, part of Springer Nature 2019

Abstract

Background Cervical laminectomy is a common strategy to decompress the spinal canal.

Methods The anatomy of the cervical spine and surrounding critical structures as viewed from the posterior approach is described. The use of a high-speed drill with a footplate attachment to make laminar troughs with an en-bloc subaxial cervical laminectomy is described with a discussion on surgical technique and complication avoidance.

Conclusion This technique allows for a safe, comfortable, and rapid decompression of the cervical spine with minimal risk. For routine cases, this may potentially be more safe and cost-effective than using a cutting bur or bone scalpel attachment.

Keywords Cervical spondylotic myelopathy · Cervical laminectomy · Adult craniotome with footplate attachment · Decompression · Durotomy · High-speed drill

Abbreviations

IAP Inferior articulating process

Relevant surgical anatomy

The cervical lamina and spinous process form the posterior bony elements of the cervical spinal canal [10]. Lamina slopes downward toward the facet and lateral masses. The lamina-facet junction is marked by a groove that identifies the lateral border of the spinal canal. The inferior articulating process (IAP) is often seen as the exposure is carried further laterally.

The posterior approach is a familiar exposure to the cervical spine for surgeons and is commonly employed when the pathology spans multiple levels. Caution must be maintained at the interlaminar spaces to avoid neural injury or a durotomy (Fig. 1).

Description of the technique

Exposure

Finding and exploiting the midline avascular plane of the ligamentum nuchae on initial dissection facilitate exposure with optimal hemostasis. Going lateral to this risks entering the paraspinal musculature with resulting bleeding. Subperiosteal exposure is continued to identify the lateral border of the lateral mass. All soft tissue should be removed from the bony elements to facilitate drilling the laminar troughs.

Surgical instrument—tapered drill bit with footplate attachment

Some surgical instruments that have been previously described include rongeurs, high-speed drill with a cutting bur attachment and an ultrasonic bone scalpel [1, 2, 6–9]. The adult-sized craniotome is comprised of a footplate attachment with a tapered router drill bit. We prefer this in adults, as

This article is part of the Topical Collection on *Spine - Other*

Electronic supplementary material The online version of this article (<https://doi.org/10.1007/s00701-019-04158-y>) contains supplementary material, which is available to authorized users.


✉ Kyle B. Mueller
kyle.mueller85@gmail.com

¹ Georgetown University Medical Center, Department of Neurosurgery, Washington, DC, USA

² Medstar Musculoskeletal Research Center, Department of Orthopaedic Surgery, Medstar Union Memorial Hospital, Baltimore, MD, USA

³ Medstar Georgetown University Hospital, Department of Surgical Services, Washington, DC, USA

Fig. 1 Key anatomical landmarks. **a** Posterior view of the bony elements of the cervical spine. The lateral masses form a box. **b** Magnified view of the lamina (star) sloping downward toward the lamina-facet junction (dashed line). The inferior articulating process is usually seen as part of the facet joint (circle). The lateral mass (asterisk) form a box with the lamina-facet junction being the medial border



compared with the pediatric craniotome, due to the more robust bone anatomy [3]. This attachment is standard in a variety of sets and has several lengths. It can be re-sterilized for multiple uses. At 2.3 mm, its size is comparable to a 2-mm Kerrison rongeur. The footplate attachment protects the soft tissue and underlying neural elements when cutting through bone. Furthermore, the router does the cutting of the bone and minimizes local heating to the surrounding neural structures and limits irrigation to preserve a clear operative field (Fig. 2).

Initial Laminotomy


It is standard practice at our institution to stabilize with instrumentation following a multi-level posterior cervical laminectomy. As such, the pilot holes for the lateral mass screws are initially drilled with a 3-mm matchstick and the lateral mass screws are either placed at this point or following the decompression. The 3-mm matchstick is then used to create a laminotomy on each side of the caudal level. The laminotomy is drilled at the caudal end of the lamina because

the ligamentum flavum, lying underneath, serves to minimize neural injury. The ligamentum flavum attaches to the middle of the lamina, and, as such, the resulting rostral aspect of the lamina has minimal soft tissue (mainly epidural fat) between the bony elements and the neural structures. Kerrison rongeurs can be further used to widen the laminotomy to facilitate placement of the footplate attachment. Curettes should be used to dissect the ligamentum flavum from the underlying bone. Preservation of the ligamentum flavum prior to drilling the laminar toughs is key to protect the underlying neural structures (Fig. 3).

Bilateral laminar troughs and en-bloc laminectomy

The footplate is placed between the lamina and ligamentum flavum. The attachment allows for better surgeon comfort than with traditional cutting tips or an ultrasonic bone scalpel. The footplate is held orthogonal to the slope of the lamina. The cut line must not extend beyond the beginning of the lateral mass so as to avoid injury to the exiting nerve roots (Fig. 4). The

Fig. 2 Surgical instrument—tapered drill bit with footplate attachment. The footplate, at 2.3 mm, is slightly larger than a 2-mm Kerrison rongeur. The drill bit used is tapered






Fig. 3 Initial laminotomy. The initial laminotomy is performed bilaterally with a 3-mm matchstick. **a** Before and **b** after completion of a laminotomy. The ligamentum flavum is under the inferior half of the

lamina and a small laminotomy is created. A 2-mm Kerrison rongeur can be used to widen the laminotomy to allow for the insertion of the footed attachment

attachment is held comfortably and the tendency to pull up with the attachment must be avoided. Doing this minimizes the ability of the drill to cut and progress forward simultaneously. The footplate will automatically move to the subsequent levels. The cutting nature of the router also allows for the added benefit of hemostasis while the cuts are being made. Furthermore, because rapid and efficient laminar troughs can be made, the use of irrigation is able to be minimized and the clear operative field is preserved. Separating the interspinous ligaments at the rostral and caudal extent of decompression then allows for an en-bloc laminectomy (Fig. 5). When removing the lamina en-bloc, be careful not to toggle or allow the bone fragments to recoil back and cause neural injury. At our institution, this bone is further used as autograft to facilitate arthrodesis.

Indications

A cervical laminectomy is useful to decompress multiple segments of the cervical spine for degenerative conditions when an anterior approach is not indicated. Moreover, its use is also advantageous for intradural tumors, vascular malformations, and continuous ossification of the posterior longitudinal ligament. The technique is also valuable for laminoplasty in certain indications.

Limitations

This technique describes an efficient method of creating laminar troughs to facilitate a multi-level posterior





Fig. 4 Bilateral laminar troughs. **a** The footplate is placed orthogonal to the lamina. The cut line should not deviate beyond the start of the lateral mass. **b** The en-bloc laminectomy is lifted to demonstrate how the

footplate fits underneath the lamina, with the footplate protecting the underlying soft tissue and thecal sac

Fig. 5 En-bloc laminectomy. **a** Bilateral troughs have been created and the interspinous ligaments cut, all the lamina can be removed en-bloc. **b** The en-bloc laminectomy has been performed. The thecal sac (asterisk) is well decompressed, and the lateral mass is preserved (box). **c** The completed lamina is removed en-bloc and can be used for auto-graft or used later for a laminoplasty



cervical laminectomy for degenerative cases. The use of more advanced osteotomy cuts for cervical deformity requires other instruments. Cervical laminectomy via a posterior approach results in limited access to anterior pathology. An increase in neck pain secondary to denervation of posterior musculature can be seen as a result of the procedure.

How to avoid complications

The initial laminotomy should be sufficient to allow for the insertion of the footplate to slide under the remaining bony lamina. Maintaining this position and allowing the footplate to “do the work” as opposed to trying to force it will allow rapid cutting of laminar troughs. Avoid the tendency to pull up as this can limit forward advancement.

Specific perioperative considerations

If there is any loss of intraoperative neuromonitoring signals along with postoperative neurological deficits, the patients’ blood pressure should be augmented to maintain perfusion postoperatively with careful monitoring of the patient’s neurological status. Any incidental durotomy should be addressed prior to closure [4, 5]. A cervical collar is often applied following the procedure. For patients at high-risk for wound healing, prophylactic negative pressure therapy can be used to reduce the risk for developing a surgical site infection.

Specific information to give to the patient about surgery and potential risks

The need for decompressing the spinal cord must be explained to the patient. If there is preoperative neurological deficit, it should

be communicated clearly that the primary purpose of surgery is to prevent further loss of function and damage to the spinal cord.

Summary of Key Points

1. Dissection along the avascular ligamentum nuchae facilitates exposure of the posterior cervical spine while minimizing blood loss.
2. The cervical lamina slopes downward toward the facet/lateral mass complex. Identifying the lamina-facet junction is a key marker for the lateral boundaries of the spinal canal.
3. Complete exposure requires meticulous dissection of the entire posterior bony elements.
4. Lateral mass screws are placed either before or after decompression (surgeon preference).
5. Initial laminotomy should be sufficient to facilitate footplate placement between the bone and ligamentum flavum.
6. The ligamentum flavum should be dissected from the underlying bone, but not removed so as to protect the underlying neural structures.
7. The footplate should be aligned orthogonal to the bony lamina at all times.
8. The footplate cut lines should not extend beyond the beginning of the lateral mass so as to avoid injury to the exiting nerve roots.
9. The footplate attachment protects the soft tissue and underlying neural elements when cutting through bone.
10. Use of the footplate attachment with a tapered router improves surgeon comfort and may be useful in resident training.

Acknowledgments We applied the sequence determines credit approach for the sequence of authors.

Compliance with ethical standards

Conflict of interest The authors declare that they have no conflict of interest.

References

1. Bydon M, Xu R, Papademetriou K, Sciubba DM, Wolinsky JP, Witham TF et al (2013) Safety of spinal decompression using an ultrasonic bone curette compared with a high-speed drill: outcomes in 337 patients. *J Neurosurg Spine* 18(6):627–633
2. Bydon M, Macki M, Xu R, Ain MC, Ahn ES, Jallo GI (2014) Spinal decompression in achondroplastic patients using high-speed drill versus ultrasonic bone curette: technical note and outcomes in 30 patients. *J Pediatr Orthop* 34(8):780–786
3. Elia C, Hariri OR, Duong J, Dong F, Sweiss R, Miulli D (2018) Use of a pediatric craniotome drill for cervical and thoracic spine decompression: safety and efficacy. *World Neurosurg* 113:e486–e489
4. Fehlings MG, Smith JS, Kopjar B, Arnold PM, Yoon ST, Vaccaro AR et al (2012) Perioperative and delayed complications associated with the surgical treatment of cervical spondylotic myelopathy based on 302 patients from the AOSpine North America Cervical Spondylotic Myelopathy Study. *J Neurosurg Spine* 16(5):425–432
5. Guerin P, El Fegoun AB, Obeid I, Gille O, Lelong L, Luc S et al (2012) Incidental durotomy during spine surgery: incidence, management and complications. A retrospective review. *Injury* 43(4):397–401
6. Ito K, Ishizaka S, Sasaki T, Miyahara T, Horiuchi T, Sakai K et al (2009) Safe and minimally invasive laminoplastic laminotomy using an ultrasonic bone curette for spinal surgery: technical note. *Surg Neurol* 72(5):470–475
7. Nakagawa H, Kim SD, Mizuno J, Ohara Y, Ito K (2005) Technical advantages of an ultrasonic bone curette in spinal surgery. *J Neurosurg Spine* 2:431–435
8. Onen MR, Yuvruk E, Akay S, Naderi S (2015) The reliability of the ultrasonic bone scalpel in cervical spondylotic myelopathy: a comparative study of 46 patients. *World Neurosurg* 84(6):1962–1967
9. Takenaka S, Hosono N, Mukai Y, Miwa T, Fuji T (2013) The use of cooled saline during bone drilling to reduce the incidence of upper-limb palsy after cervical laminoplasty. *J Neurosurg Spine* 19(4):420–427
10. Tan LA, Riew KD, Traynelis VC (2017) Cervical spine deformity-part 3: posterior techniques, clinical outcome, and complications. *Neurosurgery* 81(6):893–898

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Terms and Conditions

Springer Nature journal content, brought to you courtesy of Springer Nature Customer Service Center GmbH (“Springer Nature”).

Springer Nature supports a reasonable amount of sharing of research papers by authors, subscribers and authorised users (“Users”), for small-scale personal, non-commercial use provided that all copyright, trade and service marks and other proprietary notices are maintained. By accessing, sharing, receiving or otherwise using the Springer Nature journal content you agree to these terms of use (“Terms”). For these purposes, Springer Nature considers academic use (by researchers and students) to be non-commercial.

These Terms are supplementary and will apply in addition to any applicable website terms and conditions, a relevant site licence or a personal subscription. These Terms will prevail over any conflict or ambiguity with regards to the relevant terms, a site licence or a personal subscription (to the extent of the conflict or ambiguity only). For Creative Commons-licensed articles, the terms of the Creative Commons license used will apply.

We collect and use personal data to provide access to the Springer Nature journal content. We may also use these personal data internally within ResearchGate and Springer Nature and as agreed share it, in an anonymised way, for purposes of tracking, analysis and reporting. We will not otherwise disclose your personal data outside the ResearchGate or the Springer Nature group of companies unless we have your permission as detailed in the Privacy Policy.

While Users may use the Springer Nature journal content for small scale, personal non-commercial use, it is important to note that Users may not:

1. use such content for the purpose of providing other users with access on a regular or large scale basis or as a means to circumvent access control;
2. use such content where to do so would be considered a criminal or statutory offence in any jurisdiction, or gives rise to civil liability, or is otherwise unlawful;
3. falsely or misleadingly imply or suggest endorsement, approval, sponsorship, or association unless explicitly agreed to by Springer Nature in writing;
4. use bots or other automated methods to access the content or redirect messages
5. override any security feature or exclusionary protocol; or
6. share the content in order to create substitute for Springer Nature products or services or a systematic database of Springer Nature journal content.

In line with the restriction against commercial use, Springer Nature does not permit the creation of a product or service that creates revenue, royalties, rent or income from our content or its inclusion as part of a paid for service or for other commercial gain. Springer Nature journal content cannot be used for inter-library loans and librarians may not upload Springer Nature journal content on a large scale into their, or any other, institutional repository.

These terms of use are reviewed regularly and may be amended at any time. Springer Nature is not obligated to publish any information or content on this website and may remove it or features or functionality at our sole discretion, at any time with or without notice. Springer Nature may revoke this licence to you at any time and remove access to any copies of the Springer Nature journal content which have been saved.

To the fullest extent permitted by law, Springer Nature makes no warranties, representations or guarantees to Users, either express or implied with respect to the Springer nature journal content and all parties disclaim and waive any implied warranties or warranties imposed by law, including merchantability or fitness for any particular purpose.

Please note that these rights do not automatically extend to content, data or other material published by Springer Nature that may be licensed from third parties.

If you would like to use or distribute our Springer Nature journal content to a wider audience or on a regular basis or in any other manner not expressly permitted by these Terms, please contact Springer Nature at

onlineservice@springernature.com